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Class I medical device, manufactured by Contel S.r.l. at its headquarters Strada delle Brigole, 46 23877 Paderno d'Adda (LC) Italy in accordance with EU Directive 745/2017 MDR according to which these instructions are drawn up and are intended to ensure an adequate and safe use of the product. Contel S.r.l. is a company with ISO 9001 certified quality system

USER GUIDE FOR WRIST BRACE 1° RAY WITH TOTAL CONTROL

PK09

Thank you for choosing NewellMed®. Read the following instructions carefully for proper use of the device

THERAPEUTIC INDICATIONS

Recommended in post-surgical period of carpal tunnel and subluxations of the first metacarpo-phalangeal joint

PRECAUTIONS FOR USE

- Place the support without exerting excessive pressure and avoiding the application on parts that present injuries or excessive swelling. Do not over tighten the support so as not to cause local pressures and compressions to blood vessels and/or nerves
- For any application doubt contact your physiotherapist/ doctor/ orthopedic technician
- We recommend the replacement of the device when there is wear or even partial wear
- Wash the product regularly to ensure it remains in a hygienic condition

WARNINGS

The support is recommended if prescribed by a doctor or a physiotherapist; it must be applied by a physiotherapist/ doctor/ orthopedic technician with the utmost care to ensure proper use. Do not alter the adjustments made by the physiotherapist/ doctor/ orthopedic technician. The company declines any responsibility if the subject makes an improper use of the support or if it is used by several people. The support is not transferable as required by MDR regulation. In case of malfunction or the appearance of pain/swelling/stasis, contact the orthopaedic physiotherapist/doctor/technician. In case of particular seriousness, appropriate reporting to the manufacturer and the competent state bodies

APPLICATION

1. Open the breakaway tip
2. Insert the thumb into the thumb sleeve
3. Close the brace at the wrist
4. Adjust the thumb fastener
5. Stabilise control of the device with the lap strap
6. Ability to vary the angle of the first beam via the insert between thumb and forefinger
7. To remove the brace open the straps and pull out

MAINTENANCE

-  Wash at 30°C using neutral soap
-  Do not bleach
-  Do not tumble dry and keep away from heat sources
-  Do not iron
-  Do not dry clean

DISPOSAL

Non-biodegradable device, dispose together with normal municipal solid waste. Check the provisions of your municipality

Class I medical device, manufactured by Contel S.r.l. at its headquarters Strada delle Brigole, 46 23877 Paderno d'Adda (LC) Italy in accordance with EU Directive 745/2017 MDR according to which these instructions are drawn up and are intended to ensure an adequate and safe use of the product. Contel S.r.l. is a company with ISO 9001 certified quality system

USER GUIDE FOR WRIST BRACE 1° RAY WITH TOTAL CONTROL

PK09

Thank you for choosing NewellMed®. Read the following instructions carefully for proper use of the device

THERAPEUTIC INDICATIONS

Recommended in post-surgical period of carpal tunnel and subluxations of the first metacarpo-phalangeal joint

PRECAUTIONS FOR USE

- Place the support without exerting excessive pressure and avoiding the application on parts that present injuries or excessive swelling. Do not over tighten the support so as not to cause local pressures and compressions to blood vessels and/or nerves
- For any application doubt contact your physiotherapist/ doctor/ orthopedic technician
- We recommend the replacement of the device when there is wear or even partial wear
- Wash the product regularly to ensure it remains in a hygienic condition

WARNINGS

The support is recommended if prescribed by a doctor or a physiotherapist; it must be applied by a physiotherapist/ doctor/ orthopedic technician with the utmost care to ensure proper use. Do not alter the adjustments made by the physiotherapist/ doctor/ orthopedic technician. The company declines any responsibility if the subject makes an improper use of the support or if it is used by several people. The support is not transferable as required by MDR regulation. In case of malfunction or the appearance of pain/swelling/stasis, contact the orthopaedic physiotherapist/doctor/technician. In case of particular seriousness, appropriate reporting to the manufacturer and the competent state bodies

APPLICATION

1. Open the breakaway tip
2. Insert the thumb into the thumb sleeve
3. Close the brace at the wrist
4. Adjust the thumb fastener
5. Stabilise control of the device with the lap strap
6. Ability to vary the angle of the first beam via the insert between thumb and forefinger
7. To remove the brace open the straps and pull out

MAINTENANCE

-  Wash at 30°C using neutral soap
-  Do not bleach
-  Do not tumble dry and keep away from heat sources
-  Do not iron
-  Do not dry clean

DISPOSAL

Non-biodegradable device, dispose together with normal municipal solid waste. Check the provisions of your municipality